THUASNE

Genu Ligaflex® ROM

fr	Genouillère ligamentaire articulée avec contrôle de flexion/extension	6
en	Hinged ligament knee brace with flexion/extension control.....	8
de	Knieführungssorthese mit Gelenk und Beuge- und Streckbegrenzung	11
nl	Scharnierende ligamentaire kniebrace met buig-/strekcontrole.....	14
it	Ginocchiera per legamenti articolata con controllo di flessione/estensione.....	17
es	Rodillera ligamentaria articulada con control de flexión/extensión.....	20
pt	Joelheira ligamentar articulada com controlo de flexão/extensão	23
da	Artikuleret knæbind til ledbånd med kontrol af fleksion/ekstension.....	26
fi	Nivelletty polven nivelsidetuki, jossa on fleksion ja ekstension rajoitusmahdollisuus.....	29
sv	Ledat knästöd för ligamentskada med flexions-/extensionskontroll	31
el	Επιγονατίδα συνδέσμων αρθρωτή με έλεγχο της κάμψης/έκτασης.....	34
cs	Kolenní vazová ortéza s kloubem, umožňující nastavení flexe/extenze.....	37
pl	Przegubowa orteza stawu kolanowego z kontrolą zgięcia/wyprostu.....	40
lv	Ceļa locītavas saišu ortoze ar eņģēm un saliekšanas/atliekšanas kontroli.....	43
lt	Lankstus raiščių antkelis su sulenkimo / ištiesimo kontrole.....	46
et	Painutust-sirutust reguleerida võimaldav liigendiga põlvvetugi.....	49
sl	Pregibna opornica za koleno z nadzorom upogibanja in iztegovanja.....	51
sk	Kolenná ortéza na podporu väzov s nastavením flexie/extenzie.....	54
hu	Ízületes térdrögzítő hajlítás-nyújtás szabályozással.....	57
bg	Съчленена на коленка за връзките на коляното с контрол на сгъване/разгъване	60
ro	Orteză de genunchi ligamentară cu articulație și control al flexiei/extensiei.....	63
ru	Лигаментарный ортез для коленного сустава с шарнирами и регуляцией угла сгибания и разгибания	66
hr	Zglobna ortoza za ligamente koljena s kontrolom fleksije/ekstenzije	69
zh	带屈曲/伸展控制关节的韧带护膝.....	71
ar	داعم مفصلي للركبة وأربطتها مع ميزة التحكم بالثبي والمد.....	76

	1	32 - 34 cm
	2	35 - 37 cm
	3	38 - 41 cm
	4	42 - 45 cm
	5	46 - 49 cm
	6	50 - 53 cm


fr	Circonférence du genou
en	Knee circumference
de	Knieumfang
nl	Omvang van de knie
it	Circonferenza del ginocchio
es	Circunferencia de la rodilla
pt	Circunferência do joelho
da	Omkreds af knæ
fi	Polven ympärys
sv	Omkrets runt knä
el	Περίμετρος του γόνατος
cs	Obvod kolena
pl	Obwód kolana
lv	Celgala apkārtmērs
lt	Kelio apimtis
et	Pölvie ümbermõõt
sl	Obseg kolena
sk	Obvod kolena
hu	Térd körfogata
bg	Обиколка на коляното
ro	Circumferința genunchiului
ru	Диаметр окружности колена
hr	Opseg koljena
zh	膝关节周长
ar	محيط الركبة

fr	Stabilisation	Contrôle du mouvement	Effet proprioceptif
en	Stabilisation	Motion control	Proprioceptive effect
de	Stabilisierung	Bewegungsführung	Propriozeptive Wirkung
nl	Stabilisatie	Bewegingscontrole	Proprioceptieve werking
it	Stabilizzazione	Controllo del movimento	Effetto proprioiettivo
es	Estabilización	Control del movimiento	Efecto propioceptivo
pt	Estabilização	Controlo do movimento	Efeito propioceptivo
da	Stabilisering	Kontrol af bevægelsen	Proprioceptisk effekt
fi	Stabiloi	Liikkeen hallinta	Liikeaistiin kohdistuva vaikutus
sv	Stabilisering	Rörelsekontroll	Proprioceptiv effekt
el	Σταθεροποίηση	Έλεγχος της κίνησης	Αποτέλεσμα ιδιοδεκτικότητας
cs	Stabilizace	Kontrola pohybu	Proprioceptivní účinek
pl	Stabilizacja	Kontrola ruchu	Propriocepcja
lv	Stabilizēšana	Kustības kontrole	Proprioceptīvs efekts
lt	Stabilizavimas	Judėsio kontrolė	Propriocepčinis poveikis
et	Stabiliseerimine	Liikuvuse kontroll	Proprioseptiivne toime
sl	Stabiliziranje	Nadzor gibanja	Proprioceptivni učinek
sk	Stabilizácia	Vymedzenie rozsahu pohybu	Proprioceptívny účinok
hu	Stabilizálás	Mozgásszabályozás	Proprioceptív hatás
bg	Стабилизация	Контрол на движението	Проприоцептивен ефект
ro	Stabilizare	Controlul mișcării	Efect proprioceptiv
ru	Стабилизация	Контроль движения	Проприоцептивный эффект
hr	Stabilizacija	Kontrola pokreta	Učinak proprioceptivnog vježbanja
zh	稳定	运动控制	本体感
ar	الاستقرار	التحكم في الحركة	مفعول استقبال الحس العميق

1

MD

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en HINGED LIGAMENT KNEE BRACE WITH FLEXION/ EXTENSION CONTROL

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Bilateral product. Available in 6 sizes.

The device is composed of (see Fig. 1):

- an anatomically-shaped knitting integrating a light and soft knit at the back of the knee,
- 2 rigid hinged uprights providing knee stability ❶,
- the TM5+ hinge mimicking the knee's natural movement ❷, protected by a cover,
- 4 anterior half-straps ❸, ❹, ❺, ❻ and 4 posterior half-straps Ⓜ, Ⓨ, Ⓩ, ⓓ, ⓔ, ⓕ, ⓖ (including a suspension strap Ⓢ),
- clip buckle system ❶.

The flexion/extension adjustment system is designed to be used without tools. Extension and flexion stops are included in a plastic bag in the brace's box.

Composition

Textile components: polyamide - polyester - polyurethane - elastane.
Rigid components: polyamide - silicone.

Properties/Mode of action

Stabilization of knee joint ligaments provided by the rigid hinged side uprights.

- Extension can be adjusted to 0°, 5°, 10°, 15°, 20°, 30° and 40°.
- Flexion can be adjusted to 0°, 30°, 45°, 60°, 75° and 90°.

The brace is maintained on the leg thanks to:

- the Suspension Strap Ⓢ,
- the silicone-coated threads at the top of the knee brace.

Model totally open, making it easy to fit.

Indications

Conservative treatment of knee ligament injuries and/or ruptures (cruciate and/or lateral ligaments).

Post-operative rehabilitation.

Joint instability/laxity.

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

History of venous or lymphatic disorders.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application. Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient. It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional. Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Do not wear the product in a medical imaging device.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

To avoid damage to the knit fabric, please make sure that the self-gripping pads are always fastened to the designated self-gripping surfaces when putting it on and taking it off.

- To position the knee brace properly, place the fine knit ❶ at the back of the knee, with your leg straight; hold the uprights on either side at the hinge and perform a flexion/extension movement.

- Cover the leg with the knee brace by fastening the touch fasteners at the top and bottom.

The straps are numbered using a perforation system.

Tighten the straps following the order of the perforations, starting by the posterior strap:

- first tighten and fasten the suspension strap ❷ located above the calf,

- clip and/or tighten straps ❸, ❹, ❺, ❻, ❼, ❽ and ❾ successively.

Tightening the front/back straps adjusts the position of the uprights on the leg to optimize comfort.

If you feel the need to tighten the straps during an activity, make sure you tighten them following the steps indicated in the fitting instructions.

To make the knee brace more comfortable and improve its hold, the uprights can be shaped.

Apply a slight pressure in the desired direction.

Repeat the process on the other upright if needed.

If subsequent shaping of the uprights is needed, consult a healthcare professional.

Flexion/extension control

The adjustment of the flexion/extension must be decided and performed by the healthcare professional, not by the patient.

The product's default position is limitation of extension to 0°.

To change this setting, please follow the following instructions, to be repeated with the same limitation on both hinges.

The extension and flexion limitations are situated on plastic pieces contained in the box (Fig. 2-3).

Extension can be limited to 0°, 5°, 10°, 15°, 20°, 30° and 40°.

Flexion can be limited to 0°, 30°, 45°, 60°, 75° and 90°.

Adjusting the extension limitation:

- Choose the desired extension limitation on the corresponding plastic piece (Fig. 2).

- Open the hinge cap by using the latch located at the back; push the latch inwards and then lift the cap.

- Flex the hinge slightly and remove the existing extension limitation.

Take care to note the direction this limitation is facing.

- Insert the new extension limitation.

- Take care to check that the shape of this part follows the contours of the cap (Fig. 2).

Straighten the hinge to full extension to confirm the extension limitation is correctly positioned.

- Close the cap, and then do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

Adjusting the flexion limitation:

- Choose the desired flexion limitation on the corresponding plastic piece (Fig. 3).

- Open the hinge cap by using the latch located at the back; push the latch inwards and then lift the cap.

- Put the hinge into the maximum extension position and insert the desired flexion limitation (the flexion limitation are placed at the back of the hinge).

Make sure the hole located in flexion limitation coincides with the hole situated in the hinge (metal portion) (Fig. 3).

- Close the cap, and then do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

Final locking of the hinge cap:

- Detach the plastic piece ❶ for locking the cap from its packaging (Fig. 3).

- Once the cap is closed, point the locking piece in such a manner that it fits into the notch in the cap.

- Then push the locking piece until it comes into abutment.

Care/Maintenance

Close the self-fastening tabs before washing.

Remove lateral uprights before washing.

To do this:

- Remove the straps, making a note of their original positions.

- Pull on the top of the product to disengage the upper end of the lateral uprights from of their casings.

- Remove the uprights via the top opening.

Replace after washing:

- Insert the hinge via the top opening, with the domed part outwards and flexed towards the back.

- Insert the upper and lower ends of the uprights into their casings.

- Check that the loops protrude properly through the openings.

- Replace the straps, according to their position shown (Fig. 1).

Machine washable at 30 °C (delicate programme). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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- napsautaa ja/ tai kiristä sen jälkeen peräkkäin hihnat **©**, **Ⓞ**, **Ⓢ**, **Ⓣ**, **Ⓤ** ja **Ⓥ**.